

Product Information

Rat Bone Marrow Derived Macrophage (rBMDM)


Catalog Number	10RA-039	Cell Number	5 million cells/vial
Species	<i>Rattus norvegicus</i>	Storage Temperature	Liquid Nitrogen

Description

Macrophages are a type of white blood cell differentiated from circulating bone marrow-derived monocytes. Macrophages are responsible for detecting, engulfing and digesting cellular debris, apoptotic cells, and invading pathogens in a process called phagocytosis. Macrophages can be identified by several specific cell surface proteins including CD11b, CD11c, CD14, F4/80 (mice)/EMR1 (human), MAC-1/MAC-3, and CD68 by immunohistochemistry or flow cytometry analysis [1]. Bone marrow-derived macrophages are suitable for numerous applications including phagocytosis, gene expression profiling, and Morphological examination of cytopins using histological stains (e.g., May-Grünwald-Giemsa staining), etc.[2].

iXCells Biotechnologies provides Rat Bone Marrow Derived Macrophage (rBMDM), which were differentiated in the presence of M-CSF using the bone marrow cells isolated from adult Sprague Dawley Rat (Figure 1). rBMDM are harvested at P0 and delivered freshly or frozen. For frozen cells, each vial contains ≥ 5 million cells in 1 mL volume. rBMDM are characterized by flow cytometry analysis with $>95\%$ CD11b/CD11c+ purity (Figure 1). rBMDM are negative for mycoplasma, bacteria, yeast, and fungi. It is recommended to use Macrophage Culture Medium (Cat # MD-0097) for *in vitro* culturing of rBMDM. However, rBMDM are not recommended for long-term cultures or further expansion.

A



B




Figure 1. Rat Bone Marrow Derived Macrophage (rBMDM). (A) Phase contrast images. (B) The purity of the rBMDM was measured using CD11b/CD11c+ antibody by flow cytometry analysis.

Product Details

Tissue	Rat bone marrow
Package Size	5 million cells/vial
Passage Number	P0
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Macrophage Culture Medium (Cat# MD-0097)

Protocols

Thawing of Frozen Cells

1. Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for 1-2 minutes. Keep the cap out of water to minimize the risk of contamination.
3. Pipette the cells into a 15 mL conical tube with 5 mL fresh **Macrophage Culture Medium** (Cat# MD-0097).
4. Centrifuge at ~350 *g* for 5 minutes under room temperature.
5. Remove the supernatant and resuspend the cells in fresh Macrophage Culture Medium.
6. Culture the cell in tissue culture T75 flask or the desired culture vessel.

Safety Precaution: *it is highly recommended that protective gloves and clothing should be used when handling frozen vials.*

Standard Culture Procedure

1. rBMDM can be subcultured in **Macrophage Culture Medium** (Cat# MD-0097).
2. When cells reach ~90% confluence, remove the medium, and wash once with sterile PBS (5 mL/T75 flask).
3. Add ~3-5 mL of Cellstripper (Corning, Cat# 25-056-CI) into the flask and incubate for ~5 minutes at 37°C.

4. Add 5 mL Macrophage Culture Medium (Cat# MD-0097) into the flask and pipetting gently to detach the remaining attached cells. Collect all the cells and transfer them into a 15 mL tube.
5. Centrifuge at ~350g for 5min and resuspend the cells in desired volume of medium.
6. Seed the cells in the new culture vessels at $0.3-3 \times 10^5$ cells/cm² or other desired densities.

References

- [1]. Cline MJ and Sumner MA. Bone Marrow Macrophage Precursors. I. Some Functional Characteristics of the Early Cells of the Mouse Macrophage Series. *Blood*. 1972; 40: 62-69.
- [2]. Wynn TA and Barron L. Macrophages: master regulators of inflammation and fibrosis. *Semin Liver Dis*. 2010; 30: 245-57.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures.
© iXCells Biotechnologies 2015. All rights reserved.