

Cortical Neuron Maintenance Medium

Cortical Neuron Maintenance Medium is designed for in vitro culture of iPSC-derived cortical neurons. This medium is fully supplemented. This medium is sterile filtered and is ready for use.

| | |
|------------------------|--|
| Catalog No: | MD-0093 |
| Component: | Cortical Neuron Maintenance Medium (MD-0093): 100 mL |
| Quantity: | 100 mL |
| Product Use: | Each lot of Cortical Neuron Maintenance Medium is designed for in vitro culture of iPSC-derived cortical neurons. This medium is fully supplemented. This medium is sterile filtered and is ready for use. |
| Sterility Test: | Negative for bacteria, yeast, fungi, and mycoplasma. |
| Shipping: | Dry ice |
| Storage: | Upon delivery, store Cortical Neuron Maintenance Medium at -20°C. Protect from light. The kit is suitable for 12 months at -20°C, or 6 months at 4°C. |

Disclaimers

This product is for laboratory research purposes only and should not be used in humans. While iXCells Biotechnologies strives to provide accurate and current information on this product sheet, we do not guarantee its accuracy. Citations from scientific literature and patents are for informational purposes and may not be confirmed as accurate.

You are responsible for the safe storage, handling, and usage of this product, and iXCells Biotechnologies is not liable for any damages or injuries resulting from its receipt or use. While we make reasonable efforts to ensure the authenticity and reliability of deposited strains, we are not liable for damages arising from misidentification or misrepresentation of cultures.

Copyright

© iXCells Biotechnologies USA, Inc. 2023. All rights reserved.